禹重科技<sup>®</sup> ÜZONGLAB

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**Reference:** Rapid discrimination of Citrus reticulata ‘Chachi’ by headspace-gas chromatography-ion mobility spectrometry fingerprints combined with principal component analysis. **Food Research International**. 131 (2020) 108985.

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
Fig. 5. Topographic plot of GC-IMS spectra. (A) Pomace olive oil, (B) olive oil and (C) extra virgin olive oil.


Fig. 3. PCA results for mean-centered data for the three types of oil as obtained with UV-IMS. Samples: Pomace olive oil (red triangles); olive oil (green asterisks); extra virgin olive oil (blue squares). (For interpretation of the references to color in this figure legend, the reader is referred to the web version of the article.)

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**Reference:** Direct classification of olive oils by using two types of ion mobility spectrometers. *Anal. Chim. Acta* (2011). DOI: 10.1016/j.aca.2011.03.007.

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**Referenc:** Volatile-Compound Fingerprinting by Headspace-Gas-Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as a Benchtop Alternative to <sup>1</sup>H NMR Profiling for Assessment of the **Authenticity of honey**. **Analytical Chemistry** DOI: 10.1021/acs.analchem.7b03748 (2017)

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**Reference:** Ion mobility spectrometry coupled to gas chromatography: A rapid tool to assess eggs freshness. *Food Chemistry*, 271 (2019) 691-696.

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**Reference:** Characteristic volatiles fingerprints and changes of volatile compounds in fresh and dried *Tricholoma matsutake* Singer by HS-GC-IMS and HS-SPME-GC-MS. **Journal of Chromatography B**, 1099. (2018) 46-55.

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**Reference:** Content variations in compositions and volatile component in jujube fruits during the blacking process. *Food Sci Nutr.* 2019;1-9.

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6. Volatile-Compound Fingerprinting by Headspace Gas Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as a Benchtop Alternative to 1H NMR Profiling for Assessment of the Authenticity of Honey. **Anal. Chem.** DOI: 10.1021/acs.analchem.7b03748.
7. Untargeted and Targeted Discrimination of Honey Collected by *Apis cerana* and *Apis mellifera* Based on Volatiles Using HS-GC-IMS and HS-SPME-GC-MS. *Journal of Agricultural and Food Chemistry*. **J. Agric. Food Chem.** 2019, 67, 12144–12152.
8. Identification of terpenes and essential oils by means of static headspace gas chromatography-ion mobility spectrometry. **Anal Bioanal Chem** (2017) 409:6595–6603.
9. Volatile Metabolites of Goat Cheeses Determined by Ion Mobility Spectrometry Potential Applications in Quality Control. **Food Anal. Methods**. DOI 10.1007/s12161-014-0050-1.
10. Characterization of Arabica and Robusta Coffees by Ion Mobility Sum Spectrum. **Sensors** 2020, 20, 3123; doi:10.3390/s20113123.

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1. 气相离子迁移谱法在羊奶粉和驴肉鉴伪分析中的应用. **河北农业大学**. 2019.
2. 气相离子迁移谱分析市售燕盏挥发性物质成分. **食品工业科技**. 2020.
3. 一种鉴别卷烟真伪的方法. **专利公布号:CN 109655529 A**. 贵州中烟工业有限责任公司. 2017.
4. 一种检测菜籽粕中是否掺杂抗生素滤渣的方法. **专利公布号:CN 108693004 A**. 中国农业科学院农业质量标准与检测技术研究所. 2018.
5. Direct differentiation of herbal medicine for volatile components by a multicapillary column with ion mobility spectrometry method. **J. Sep. Sci.** 2015, 38, 3205–3208.
6. Target vs spectral fingerprint data analysis of Iberian ham samples for avoiding labelling fraud using headspace – gas chromatography–ion mobility spectrometry. **Food Chemistry**. 246 (2018) 65-73.
7. A green triple-locked strategy based on volatile-compound imaging, chemometrics, and markers to discriminate winter honey and sapium honey using headspace gas chromatography-ion mobility spectrometry. **Food Research International**. <https://doi.org/10.1016/j.foodres.2019.01.004>.
8. Untargeted headspace gas chromatography – Ion mobility spectrometry analysis for detection of adulterated honey. **Talanta**. 205 (2019) 120123.
9. Detection of Peanut Oil Adulteration Mixed with Rapeseed Oil Using Gas Chromatography and Gas Chromatography–Ion Mobility Spectrometry. **Food Analytical Methods**. <https://doi.org/10.1007/s12161-019-01571-y>.
10. Molecular Analysis for Determining Differences and Similarities of Perfume Dupes. **Global Ingredients & Formulations Gulde**. 2016.

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1. 基于气相色谱与离子迁移谱的车内气味等级评价方法. 专利公布号:CN 111007176 A. 冲之智能科技(天津)有限公司. 2020.
2. 基于气相离子迁移谱技术的葵花籽油精炼程度的研究. **食品科学**. 2019.
3. 一种中药甘草的质量评价方法及其应用. 专利公布号:CN 110749681 A. **山东省中医药研究院**. 2020.
4. A robustness study of calibration models for olive oil classification: Targeted and non-targeted fingerprint approaches based on GC-IMS. **Food Chemistry**. 288 (2019) 315-324.
5. Resolution-optimized headspace gas chromatography-ion mobility spectrometry (HS-GC-IMS) for non-targeted olive oil profiling. **Anal Bioanal Chem**. DOI 10.1007/s00216-017-0338-2.
6. Direct classification of olive oils by using two types of ion mobility spectrometers. **Anal. Chim. Acta** (2011), doi:10.1016/j.aca.2011.03.007.
7. Gas Chromatography-Ion Mobility Spectrometry Detection of Odor Fingerprint as Markers of Rapeseed Oil Refined Grade. **Journal of Analytical Methods in Chemistry**. <https://doi.org/10.1155/2019/3163204>.
8. Usefulness of GC-IMS for rapid quantitative analysis without sample treatment: Focus on ethanol, one of the potential classification markers of olive oils. **LWT - Food Science and Technology**. 120 (2020) 108897.
9. Use of headspace-gas chromatography-ion mobility spectrometry to detect volatile fingerprints of palm fibre oil and sludge palm oil in samples of crude palm oil. **BMC Res Notes** (2019) 12:229 <https://doi.org/10.1186/s13104-019-4263-7>.
10. Characteristic volatiles fingerprints and profiles determination in different grades of coconut oil by HS-GC-IMS and HS-SPME-GC-MS. **International Journal of Food Science and Technology**. 2020. doi:10.1111/ijfs.14664.

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6. Changes in the Volatile Components of Candied Kumquats in Different Processing Methodologies with Headspace-Gas Chromatography-Ion Mobility Spectrometry. **Molecules** 2019, 24, 3053; doi:10.3390/molecules24173053.
7. Evaluation by electronic tongue and headspace-GC-IMS analyses of the flavor compounds in dry-cured pork with different salt content. **Food Research International**.  
<https://doi.org/10.1016/j.foodres.2020.109456>.
8. Biochemical changes induced by dominant bacteria in chill-stored silver carp (*Hypophthalmichthys molitrix*) and GC-IMS identification of volatile organic compounds. **Food Microbiology**. 84 (2019) 103248.
9. Mechanism of aroma compounds changes from sea cucumber peptide powders (SCPPs) under different storage conditions. **Food Research International**. 128 (2020) 108757.
10. Changes in volatile flavor compounds of peppers during hot air drying process based on headspace-gas chromatography-ion mobility spectrometry (HS-GC-IMS). **J Sci Food Agric**. 2020; 100: 3087–3098.

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8. Ion mobility spectrometry coupled to gas chromatography: A rapid tool to assess eggs freshness. **Food Chemistry**. 271 (2019) 691-696.
9. Analysis of the volatile compounds associated with pickling of ginger using headspace gas chromatography - ion mobility spectrometry. **Flavour Fragr J**. DOI: 10.1002/ffj.3530.
10. Ion Mobility Spectrometry versus Classical Physico-chemical Analysis for Assessing the Shelf Life of Extra Virgin Olive Oil According to Container Type and Storage Conditions. **J. Agric. Food Chem**. 2015, 63, 2179–2188.

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1. 臭鳊鱼低温发酵过程中品质及挥发性物质的变化. **中国食品学报**. 2020.
2. 气相离子迁移谱联用技术评定大米霉变程度的研究应用. **中国粮油学报**. 2019.
3. Characterization of aroma profile and characteristic aromas during lychee wine fermentation. **J Food Process Preserv**. DOI: 10.1111/jfpp.14003.
4. A Study on Volatiles Metabolites Screening by HS-SPME-GC-MS and HS-GC-IMS for Discrimination and Characterization of White and Yellowed Rice. **Cereal Chemistry**. 2020; 00:1–9.. DOI: 10.1002/cche.10264.
5. Evaluation of fast volatile analysis for detection of Botrytis cinerea infections in strawberry. **Food Microbiology** 32 (2012) 406e414.
6. Effects of different probiotic combinations on the components and bioactivity of Spirulina. **J Basic Microbiol**. 2020;1–15.
7. Changes in volatile compounds of fermented minced pepper during natural and inoculated fermentation process based on headspace–gas chromatography–ion mobility spectrometry. **Food Sci Nutr**. 2020; 00:1–18.
8. Target identification of volatile metabolites to allow the differentiation of lactic acid bacteria by gas chromatography-ion mobility spectrometry. **Food Chemistry**. 220 (2017) 362–370.
9. Rapid detection of hydrogen sulfide produced by pathogenic bacteria in focused growth media using SHS-MCC-GC-IMS. **Microchemical Journal**. 140 (2018) 232–240.
10. Profiling and characterization of odorous volatile compounds from the industrial fermentation of erythromycin. **Environmental Pollution**. 255 (2019) 113130.